


|                            |            |       |
|----------------------------|------------|-------|
| CATALOG NUMBER             | SEE BELOW  |       |
| WORKING RANGE (μL)         | 50 – 250μL |       |
| MATERIAL                   | PS         |       |
|                            | MM         | INCH  |
| A, BASE LENGTH             | 127,7      | 5,03  |
| B, BASE WIDTH              | 85,7       | 3,37  |
| C, TOP LENGTH              | 123,5      | 4,86  |
| D, TOP WIDTH               | 81,3       | 3,20  |
| E, FLANGE LONG SIDE        | 2,7        | 0,11  |
| F, OVERALL HEIGHT          | 14,5       | 0,57  |
| G1, DIA TOP                | Ø6,6       | Ø0,26 |
| G2, DIA BOTTOM             | Ø6,1       | Ø0,24 |
| H, A - I LOCATION (X)      | 14,3       | 0,56  |
| J, A - I LOCATION (Y)      | 11,4       | 0,45  |
| K, WELL TO WELL SPACING    | 9          | 0,35  |
| L, WELL DEPTH              | 11         | 0,43  |
| M, BOTTOM OF WELL DISTANCE | 3,5        | 0,14  |
| N, FLANGE SHORT SIDE       | 2,7        | 0,11  |

| CATALOG NO | SURFACE  | COLOR | STERILE | WITH LID* | UNITS/ PK/CASE |
|------------|----------|-------|---------|-----------|----------------|
| 437796     | MAXISORP | WHITE | -       | -         | 5/60           |
| 437842     | POLYSORP | WHITE | -       | -         | 5/60           |
|            |          |       |         |           |                |

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**ThermoFisher**  
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280  
DK-4000 Roskilde • Denmark  
Tel +45 4631 2000 • Fax +45 4631 2042  
E-mail: info.nunc@thermofisher.com

[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT:

C96 MICROWELL PLATE, WHITE, PS

DRAWN BY/DATE:

CHANGE ORDER NO.:

VALID FROM:

APPROVED BY/DATE:

FORMA1

ITEM NUMBER: MATERIAL:

DRAWING NO.:

VERSION:

PAGE TYPE: SCALE:

METRIC

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC